



#3201

Agglomeration Tendency

EVALUATION OF AGGLOMERATED FRACTIONS OF CONDITIONED
POWDER SAMPLES

DESCRIPTION

Powdered samples can be conditioned at different relative humidities and temperatures and with different loads. The test allows evaluating the agglomerated mass fraction as a function of time and the storage conditions. Besides, the equilibrium and critical relative humidities can be determined.

AVAILABLE EQUIPMENT

- Moisture analyser OHAUS MB45
- Titulador Karl Fischer Mettler Toledo
- RH/T stability chambers

APPLICATIONS

Pharmaceuticals | Dietary supplements | Cosmetics | Phytomedicines |
Fertilizers | Cereals | powered food |

CONTACT US

Office of Technology Transfer

✉ ott@plapiqui.edu.ar

🌐 plapiqui.edu.ar/ott

☎ +54 291 4037200 - Int 217/214

📱 +54 9 291 4261644