



#5204

Dissolution apparatus (AGILENT, 708 DS)

DESCRIPTION

Automated system for dissolution testing of powders, granules, and solid pharmaceutical forms. It allows the evaluation of dissolution rate and active ingredient release under controlled conditions.

MAIN FEATURES

- Vessel: 8
- Media Volume: 400-1000mL
- Temperature Range: +5 to 55 °C
- Speed Range: 50-250 RPM
- Accessories: Basket, paddle
- Dissolution Sampling Station

MAIN APPLICATIONS

- Dissolution rate profiles of pharmaceutical formulations
- In vitro bioequivalence studies
- Pharmaceutical development.

ASSOCIATED SERVICES

3206 - 3207 - 3209

CONTACT US

Office of Technology
Transfer

✉ ott@plapiqui.edu.ar

🌐 plapiqui.edu.ar/ott

☎ +54 291 4037200 - Int 217/214

📠 +54 9 291 4261644