



#5203

Tablet disintegration tester (SCOUT)

DESCRIPTION

Equipment for determining the disintegration time of solid pharmaceutical forms such as tablets and capsules, in accordance with international pharmacopeias.

MAIN FEATURES

- 6 test stations
- Temperature: $37 \pm 0,5$ °C
- Stroke Speed: 30 ± 1 strokes/min

MAIN APPLICATIONS

- Quality control of pharmaceutical products, formulation development, and validation of formulations.

ASSOCIATED SERVICES

3204 - 3207

CONTACT US

Office of Technology
Transfer

✉ ott@plapiqui.edu.ar

🌐 plapiqui.edu.ar/ott

☎ +54 291 4037200 - Int 217/214

📠 +54 9 291 4261644